


LATEST EDITION

ORTHOPEDIC IMPLANTS



CORTICAL SCREW 3.5 mm

1.25 mm 1.75 mm PITCH ARE AVAILABLE

Same Article Available in Titanium

Article #	Length	Article #	Length
SS -700	12 mm	SS -720	32 mm
SS -702	14 mm	SS -722	34 mm
SS -704	16 mm	SS -724	36 mm
SS -706	18 mm	SS -726	38 mm
SS -708	20 mm	SS -728	40 mm
SS -710	22 mm		
SS -712	24 mm		
SS -714	26 mm		
SS -716	28 mm		
SS -718	30 mm		

CORTICAL SCREW 4.5 mm

1.25 mm 1.75 mm PITCH ARE AVAILABLE

Same Article Available in Titanium

Article #	Length	Article #	Length
SS -730	12 mm	PIO-754	36 mm
SS -732	14 mm	PIO-756	38 mm
SS -734	16 mm	PIO-758	40 mm
SS -736	18 mm	PIO-760	42 mm
SS -738	20 mm	PIO-762	44 mm
SS -740	22 mm	PIO-764	46 mm
SS -742	24 mm	PIO-766	48 mm
SS -744	26 mm	PIO-768	50 mm
SS -746	28 mm	PIO-770	55 mm
SS -748	30 mm	PIO-772	60 mm
SS -750	32 mm	PIO-774	65 mm
SS -752	34 mm	PIO-776	70 mm

MALLEOLAR SCREW 4.5 mm

Same Article Available in Titanium

Article #	Length
SS -778	25 mm
SS -780	30 mm
SS -782	35 mm
SS -784	40 mm
SS -786	45 mm
SS -788	50 mm
SS -790	55 mm
SS -792	60 mm
SS -794	65 mm
SS -796	70 mm

CANCELLOUS SCREW 4.0 mm

Same Article Available in Titanium

Article #	Length	Article #	Length
SS-798	14 mm	SS-812	28 mm
SS-800	16 mm	SS-814	30 mm
SS-802	18 mm	SS-816	35 mm
SS-804	20 mm	SS-818	40 mm
SS-806	22 mm	SS-820	45 mm
SS-808	24 mm	SS-822	50 mm
SS-810	26 mm	SS-824	60 mm



CANCELLOUS SCREW 6.5 mm X 16 mm Thread

Same Article Available in Titanium

Article #	Length	Article #	Length
SS-826	25 mm	SS-844	70 mm
SS-828	30 mm	SS-846	75 mm
SS-830	35 mm	SS-848	80 mm
SS-832	40 mm	SS-850	85 mm
SS-834	45 mm	SS-852	90 mm
SS-836	50 mm	SS-854	95 mm
SS-838	55 mm	SS-856	100 mm
SS-840	60 mm	SS-858	105 mm
SS-842	65 mm	SS-860	110 mm



CANCELLOUS SCREW 6.5 mm X 32 mm Thread

Same Article Available in Titanium

Article #	Length	Article #	Length
SS-862	45 mm	SS-876	80 mm
SS-864	50 mm	SS-878	85 mm
SS-866	55 mm	SS-880	90 mm
SS-868	60 mm	SS-882	95 mm
SS-870	65 mm	SS-884	100 mm
SS-872	70 mm	SS-886	105 mm
SS-874	75 mm	SS-888	110 mm



CANCELLOUS SCREW 6.5 mm Full Thread

Same Article Available in Titanium

Article #	Length	Article #	Length
SS -890	25 mm	SS -910	75 mm
SS -892	30 mm	SS -912	80 mm
SS -894	35 mm	SS -914	85 mm
SS -896	40 mm	SS -916	90 mm
SS -898	45 mm	SS -918	95 mm
SS -900	50 mm	SS -920	100 mm
SS -902	55 mm	SS -922	105 mm
SS -904	60 mm	SS -924	110 mm
SS -906	65 mm		
SS -908	70 mm		



LAG SCREW

Same Article Available in Titanium

Article #	Length
SS -926	50 mm
SS -928	55 mm
SS -930	60 mm
SS -932	65 mm
SS -934	70 mm
SS -936	75 mm
SS -938	80 mm
SS -940	85 mm
SS -942	90 mm
SS -944	95 mm
SS -946	100 mm
SS -948	105 mm
SS -950	110 mm



WASHERS FOR CANCELLOUS AND CORTICAL SCREW

Article # :PIO-952
Dia 40 mm
For Cancellous Screw

Article # :PIO-954
Dia 65 mm
For Cancellous Screw
& Canulated Screw



CANULATED SCREW 7.5 mm x 32 mm THREAD

Same Article Available in Titanium

Article #	Length
SS -956	55 mm
SS -958	60 mm
SS -960	65 mm
SS -962	70 mm
SS -964	75 mm
SS -966	80 mm
SS -968	85 mm
SS -970	90 mm
SS -972	95 mm
SS -974	100 mm
SS -976	105 mm
SS -978	110 mm

CANULATED SCREW 7.5 mm x 16 mm THREAD

Same Article Available in Titanium

Article #	Length
SS -980	55 mm
SS -982	60 mm
SS -984	65 mm
SS -986	70 mm
SS -988	75 mm
SS -990	80 mm
SS -992	85 mm
SS -994	90 mm
SS -996	95 mm
SS -998	100 mm
SS -1002	105 mm
SS -1004	110 mm





DHS PLATE 95°
Same Article Available in Titanium

Article #	Holes	Length
SS -1006	4 hole	70 mm
SS -1008	5hole	85 mm
SS -1010	6 hole	100 mm
SS -1012	7 hole	115 mm
SS -1014	8 hole	130 mm
SS -1016	10 hole	163 mm
SS -1018	12 hole	198 mm
SS -1020	14 hole	225 mm
SS -1022	16 hole	252 mm

JEWETT NAIL PLATE 120°
Same Article Available in Titanium

Article #	Holes	Length	Article #	Holes	Length
SS -1024	4 hole	70 mm	SS -1054	6 hole	85 mm
SS -1026	4 hole	75 mm	SS -1056	6 hole	90 mm
SS -1028	4 hole	80 mm	SS -1058	6 hole	95 mm
SS -1030	4 hole	85 mm	SS -1060	8 hole	70 mm
SS -1032	4 hole	90 mm	SS -1062	8 hole	75 mm
SS -1034	4 hole	95 mm	SS -1064	8 hole	80 mm
SS -1036	5 hole	70 mm	SS -1066	8 hole	85 mm
SS -1038	5 hole	75 mm	SS -1068	8 hole	90 mm
SS -1040	5 hole	80 mm	SS -1070	8 hole	95 mm
SS -1042	5 hole	85 mm	SS -1072	10 hole	70 mm
SS -1044	5 hole	90 mm	SS -1074	10 hole	75 mm
SS -1046	5 hole	95 mm	SS -1076	10 hole	80 mm
SS -1048	6 hole	70 mm	SS -1078	10 hole	85 mm
SS -1050	6 hole	75 mm	SS -1080	10 hole	90 mm
SS -1052	6 hole	80 mm	SS -1082	10 hole	95 mm



AUSTIN MOORE FORMAL HIP PROSTHESIS

NARROW STEM

Article #	Length	Article #	Length
SS -1084	38 mm	SS -1104	48 mm
SS -1086	39 mm	SS -1106	49 mm
SS -1088	40 mm	SS -1108	50 mm
SS -1090	41 mm	SS -1110	51 mm
SS -1092	42 mm	SS -1112	52 mm
SS -1094	43 mm	SS -1114	53 mm
SS -1096	44 mm	SS -1116	54 mm
SS -1098	45 mm	SS -1118	55 mm
SS -1100	46 mm	SS -1120	56 mm
SS -1102	47 mm	SS -1122	57 mm

REGULAR STEM

Article #	Length	Article #	Length
SS -1124	38 mm	SS -1144	48 mm
SS -1126	39 mm	SS -1146	49 mm
SS -1128	40 mm	SS -1148	50 mm
SS -1130	41 mm	SS -1150	51 mm
SS -1132	42 mm	SS -1152	52 mm
SS -1134	43 mm	SS -1154	53 mm
SS -1136	44 mm	SS -1156	54 mm
SS -1138	45 mm	SS -1158	55 mm
SS -1140	46 mm	SS -1160	56 mm
SS -1142	47 mm	SS -1162	57 mm



SEMI TUBULAR PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS -1164	3 hole	55 mm
SS -1168	4 hole	71 mm
SS -1170	5 hole	87 mm
SS -1172	6 hole	103 mm
SS -1174	7 hole	119 mm
SS -1176	8 hole	135 mm
SS -1178	9 hole	151 mm
SS -1180	10 hole	167 mm
SS -1182	12 hole	199 mm



MICRO PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS -1184	4 hole	23 mm
SS -1186	5 hole	29 mm
SS -1188	6 hole	35 mm

Long Barrel



Short Barrel



DHS PLATE 135°

Same Article Available in Titanium

Article #	Holes	Length
SS -1190	4 hole	80 mm
SS -1192	5 hole	96 mm
SS -1194	6 hole	112 mm
SS -1196	8 hole	144 mm
SS -1198	10 hole	176 mm
SS -1200	12 hole	208 mm
SS -1202	14 hole	240 mm
SS -1204	16 hole	272 mm

135°
145°
150°

Long & Short Barrel are Available in all Degrees



DCP BROAD

Same Article Available in Titanium

Article #	Holes	Length
SS-1206	5 hole	87 mm
SS-1208	6 hole	103 mm
SS-1210	7 hole	119 mm
SS-1212	8 hole	135 mm
SS-1214	10 hole	167 mm
SS-1216	12 hole	199 mm
SS-1218	14 hole	231 mm
SS-1220	16 hole	263 mm



DCP NERROW

Same Article Available in Titanium

Article #	Holes	Length
SS-1222	4 hole	71 mm
SS-1224	5 hole	83 mm
SS-1226	6 hole	103 mm
SS-1228	7 hole	119 mm
SS-1230	9 hole	135 mm
SS-1232	10 hole	167 mm
SS-1234	12 hole	199 mm
SS-1236	14 hole	231 mm



DCP SMALL FREGMENT

Same Article Available in Titanium

Article #	Holes	Length
SS-1238	5 hole	61 mm
SS-1240	6 hole	73 mm
SS-1242	7 hole	85 mm
SS-1244	8 hole	97 mm
SS-1246	10 hole	121 mm
SS-1248	12 hole	145 mm
SS-1250	14 hole	157 mm



1/3 TUBULAR PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS-1252	3 hole	37 mm
SS-1254	4 hole	49 mm
SS-1256	5 hole	61 mm
SS-1258	6 hole	73 mm
SS-1260	7 hole	85 mm
SS-1262	8 hole	97 mm
SS-1264	9 hole	109 mm
SS-1266	10 hole	121 mm



CONDYLAR BLADE PLATE 95°

Same Article Available in Titanium

Article #	Holes	Length
SS -1268	5 hole	50 mm
SS -1270	5 hole	60 mm
SS -1272	5 hole	70 mm
SS -1274	5 hole	80 mm
SS -1276	7 hole	50 mm
SS -1278	7 hole	60 mm
SS -1280	7 hole	70 mm
SS -1282		
SS -1284	9 hole	50 mm
SS -1286	9 hole	60 mm
SS -1288	9 hole	70 mm
SS -1290	9 hole	80 mm
SS -1292	12 hole	50 mm
SS -1294	12 hole	60 mm
SS -1296	12 hole	70 mm
SS -1298	12 hole	80 mm
SS -1300	14 hole	50 mm
SS -1302	14 hole	60 mm
SS -1304	14 hole	70 mm
SS -1306	14 hole	80 mm

CONDYLAR BLADE PLATE 130°

Same Article Available in Titanium

Article #	Holes	Length
SS -1308	5 hole	50 mm
SS -1310	5 hole	60 mm
SS -1312	5 hole	70 mm
SS -1314	5 hole	80 mm
SS -1316	7 hole	50 mm
SS -1318	7 hole	60 mm
SS -1320	7 hole	70 mm
SS -1322	7 hole	80 mm
SS -1324	9 hole	50 mm
SS -1326	9 hole	60 mm
SS -1328	9 hole	70 mm
SS -1330	9 hole	80 mm
SS -1332	12 hole	50 mm
SS -1334	12 hole	60 mm
SS -1336	12 hole	70 mm
SS -1338	12 hole	80 mm
SS -1340	14 hole	50 mm
SS -1342	14 hole	60 mm
SS -1344	14 hole	70 mm
SS -1346	14 hole	80 mm



CASPER PLATE & PEDICAL SCREW

Same Article Available in Titanium

PEDICAL SCREW

Same Article Available in Titanium

Article #	Length
SS -1348	30 mm
SS -1350	35 mm
SS -1352	40 mm
SS -1354	45 mm



Article #	Length
SS -1356	28 mm
SS -1358	30 mm
SS -1360	32 mm
SS -1362	37 mm
SS -1364	40 mm
SS -1366	43 mm
SS -1368	46 mm
SS -1370	48 mm
SS -1372	50 mm
SS -1374	54 mm
SS -1376	60 mm
SS -1378	66 mm
SS -1380	72 mm



SCHANZ SCREW

Same Article Available in Titanium

Article #	Size	Length
SS -1382	3.0 mm	90 to 150 mm
SS -1384	3.5 mm	90 to 150 mm
SS -1386	4.0 mm	90 to 150 mm
SS -1388	4.5 mm	90 to 150 mm
SS -1390	5.0 mm	90 to 150 mm

STEMIEN PIN

Same Article Available in Titanium

Article #	Size	Length
SS -1392	3.0 mm	90 to 150 mm
SS -1394	3.5 mm	90 to 150 mm
SS -1396	4.0 mm	90 to 150 mm
SS -1398	4.5 mm	90 to 150 mm
SS -1400	5.0 mm	90 to 150 mm

RUSH NAIL

Same Article Available in Titanium

Article #	Size	Length
SS -1402	3.0 mm	90 to 150 mm
SS -1404	3.5 mm	90 to 150 mm
SS -1406	4.0 mm	90 to 150 mm
SS -1408	4.5 mm	90 to 150 mm
SS -1410	5.0 mm	90 to 150 mm



RUSH NAIL

STEMIEN PIN

SCHANZ SCREW

Same Article Available in Titanium



SS-1412 2.0 mm



SS-1414 2.7 mm



SS-1416 3.5 mm



SS-1418 4.5 mm

RECONSTRUCTION PLATES 4.5 mm Straight

Do not bend more than 15° Fixation with 4.5 & 2.7, 3.5, 2 mm are also Available

LIMITED CONTACT DYNAMIC COMPRESSION PLATE (LC-DCP) FOR 4.5 mm SCREW (BROAD)

Fixation with 4.5 & 6.5 mm Screws. Undercuts: For reduced vascular damage. The end holes can be used for cancellous bone screws. The LC-DCP holes allow compression in both directions of the longitudinal axis. Indication: As a tension band plate on the femur and for Humeral Pseudarthrosis (not intended for the Tibia)

Thickness: 6 mm, width: 17.5 mm, Hole Spacing: 18 mm

SS-1420



SS-1422



K-NAIL

SS-1424 to SS-1478 A



SS-1424 to SS-1478 B




SS-1424 to SS-1478 C



INTERLOCKING NAIL Also Available in Titanium

SS-1424 to SS-1478 D



SS-1424 to SS-1478 E



Article #	Dia	Length	Article #	Dia	Length	Article #	Dia	Length
SS-1424	9 mm	34 cm	SS-1448	10 mm	44 cm	SS-1472	12 mm	40 cm
SS-1426	9 mm	36 cm	SS-1450	10 mm	46 cm	SS-1474	12 mm	42 cm
SS-1428	9 mm	38 cm	SS-1452	11 mm	34 cm	SS-1476	12 mm	44 cm
SS-1430	9 mm	40 cm	SS-1454	11 mm	36 cm	SS-1478	12 mm	46 cm
SS-1432	9 mm	42 cm	SS-1456	11 mm	38 cm			
SS-1434	9 mm	44 cm	SS-1458	11 mm	40 cm			
SS-1436	9 mm	46 cm	SS-1460	11 mm	42 cm			
SS-1438	10 mm	34 cm	SS-1462	11 mm	44 cm			
SS-1440	10 mm	36 cm	SS-1464	11 mm	46 cm			
SS-1442	10 mm	38 cm	SS-1466	12 mm	34 cm			
SS-1444	10 mm	40 cm	SS-1468	12 mm	36 cm			
SS-1446	10 mm	42 cm	SS-1470	12 mm	38 cm			



SMALL T PLATE OBLIQUE ANGLE 3 HOLE HEAD

Same Article Available in Titanium


Article #	Holes	Length
SS-1480	3 hole	52 mm
SS-1482	4 hole	63mm
SS-1484	5 hole	74mm



SMALL BUTTRESS PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS-1486	4 hole	54 mm
SS-1488	5 hole	66 mm
SS-1490	6 hole	78 mm
SS-1492	4 hole	54 mm
SS-1494	5 hole	66 mm
SS-1496	6 hole	78 mm



L BUTTRESS PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS-1498	4 hole	84 mm
SS-1500	5 hole	101 mm
SS-1502	6 hole	117 mm
SS-1504	4 hole	84 mm
SS-1506	5 hole	101 mm
SS-1508	6 hole	117 mm

Y PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS-1510	3 hole	77 mm
SS-1512	5 hole	106 mm
SS-1514	3 hole	77 mm
SS-1516	5 hole	106 mm





SMALL T PLATE 3 HOLE HEAD

Same Article Available in Titanium

Article #	Holes	Length
SS -1518	3 hole	44 mm
SS -1520	4 hole	56 mm
SS -1522	5 hole	68 mm
SS -1524	6 hole	80 mm



SMALL T PLATE 4 HOLE HEAD

Same Article Available in Titanium

Article #	Holes	Length
SS -1526	3 hole	44 mm
SS -1528	4 hole	56 mm
SS -1530	5 hole	68 mm
SS -1532	6 hole	80 mm



T BUTTRESS PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS -1534	3 hole	81 mm
SS -1536	4 hole	98 mm
SS -1538	5 hole	115 mm
SS -1540	6 hole	132 mm
SS -1542	7 hole	149 mm
SS -1544	8 hole	166 mm
SS -1546	10 hole	200 mm



T STANDARD PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS -1548	3 hole	81 mm
SS -1550	4 hole	98 mm
SS -1552	5 hole	115 mm
SS -1554	6 hole	132 mm
SS -1556	7 hole	149 mm
SS -1558	8 hole	166 mm
SS -1560	10 hole	200 mm



CLOVERLEAF PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS -1562	3 hole	88 mm
SS -1564	4 hole	104 mm
SS -1566	5 hole	120 mm
SS -1568	6 hole	136 mm



SPOIN PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS -1570	5 hole	100 mm
SS -1572	6 hole	120 mm



COBRA HEAD PLATE

Same Article Available in Titanium

Article #	Holes	Length
SS -1574	7 hole	157 mm
SS -1576	8 hole	170 mm
SS -1578	9 hole	186 mm
SS -1580	10 hole	202 mm
SS -1582	11 hole	218 mm
SS -1584	12 hole	234 mm

PLATES FOR DISTAL FEMUR LATERAL

For Use with 4.5 mm Screws


Condylar Buttress Plates

Primary Indication: As a buttress for femoral condyles with multiple fragments (Left and Right Plates)

Same Article Available in Titanium



Article #	Length	Holes in shaft
SS -1586	158 mm	7 Right
SS -1588	158 mm	7 Left
SS -1590	190 mm	9 Right
SS -1592	190 mm	9 Left
SS -1594	221 mm	11 Right
SS -1596	221 mm	11 Left
SS -1598	253 mm	13 Right
SS -1600	253 mm	13 Left
SS -1602	285 mm	15 Right
SS -1604	285 mm	15 Left



BONE INSTRUMENTS

SLICK SURGICO®
We care for you!



SS-1644



SS-1646
Lexer



SS-1648
Lexer



SS-1650
Lexer



SS-1652
Screw Driver



SS-1654
Screw Driver



SS-1656
Screw Driver



SS-1658
Screw Driver



SS-1660



SS-1662
Mallets



SS-1664



SS-1666
Handle with Quick
Coupling



SS-1668
Small Wire Cutter
16 cm Long



SS-1670
Wire Cutter Large
Version 22 cm Long



SS-1672
Flate Nose Parallel Plier
16 cm Long



SS-1674
Vice Grip Self Locking
18 cm Long



SS-1676
Wire Bending Plier
15.5 cm Long




SS-1678
Bone Holding Forcep



SS-1680
Reduction Forcep
with Points

SS-1682
Reduction Forcep
with Points



SS-1684
Lane Bone Holding Forcep
Ranges: 33,45 cm

SS-1686
Colclough Forceps
Ranges : 27 cm up



SS -1688




SS -1690
Hand Drill with Adjustable
Handle and 3 Jaw Chuck



SS -1694
Ralk



SS -1692
Tightening
wrench only



SS -1696
Bunnell



Please call Customer Service today for additional information, an in-office Demonstration, to place an order, or to receive other Slick Surgico catalogs.

www.slicksurgico.com | sales@slicksurgico.com